Technical standards
for
supply related samples
for NAT testing

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CONFIDENTIAL: This material is the property of Takeda and must not be disclosed except as authorized in writing.

Contents

1.0	Purpose.....	3
2.0	Approved plasma sample tubes.....	3
2.1	Sample tubes 75mm.....	3
2.2	Gel sample tubes 100mm.....	3
3.0	Barcode position on the sample tube.....	4
4.0	Packing of samples.....	4
5.0	Labels for sample tubes and sample boxes.....	5
5.1	NAT sample tubes.....	5
5.1.1	Minimum requirements of content:	5
5.2	Box for NAT sample tubes.....	5
5.2.1	Minimum requirements of content:	5
5.3	Bigger box (outer box for NAT sample tube boxes)	5
5.3.1	Minimum requirements of content:	5
6.0	Hensikt.....	6
7.0	Godkjente prøverør til plasma.....	6
7.1	Prøverør 75 mm	6
7.2	Prøverør med gel 100 mm.....	6
8.0	Plassering av strekkode på prøverør	7
9.0	Pakking av prøver	7
10.0	Etiketter til prøverør og prøveesker.....	8
10.1	NAT-prøverør.....	8
10.1.1	Minimumskrav til innhold:.....	8
10.2	Eske til NAT-prøverør.....	8
10.2.1	Minimumskrav til innhold:.....	8
10.3	Stor eske (yttereske til esker med NAT-prøverør)	8
10.3.1	Minimumskrav til innhold:.....	8

CONFIDENTIAL: This material is the property of Takeda and must not be disclosed except as authorized in writing.

1.0 Purpose

The sample is the test equivalent for the corresponding plasma donation and must therefore be clearly assignable to it.

Haemolytic, lipemic or icteric samples must not be supplied.

NAT testing sample must be made available for every plasma donation supplied to Takeda.

Samples need to be taken as part of the blood or plasma collection procedure and need to be stored according to the valid guidelines.

Any contamination or mix up between the samples of different donations must be excluded by appropriate measures.

No testing material must be removed from tubes by suppliers or external laboratories.

Storage and transport of the samples must be performed equally as the plasma donations and in compliance with the specifications of the tube manufacturer.

In case any temperature deviations occur, the contact person at Plasma Logistics & Plasma Analytics must be informed in advance in order to clarify if the samples can still be tested.

N.B.: If a sample cannot be tested, this will lead to the destruction of the corresponding plasma donation.

2.0 Approved plasma sample tubes

Black caps are not allowed.

2.1 Sample tubes 75mm

non ridged, volume 3 ml, no additives

Filling height / filling volume: min. 35 mm / corresponds to min. 2.5 ml
max. 50 mm / corresponds to max. 4 ml

- Greiner Vacuette (red cap) cat. no. 454 224 ¹⁾ or
- Greiner Vacuette (white cap) cat. no. 454 045 ¹⁾ or
- Greiner Vacuette (yellow cap) cat. no. 454 256 ¹⁾ or
- Becton Dickinson Vacutainer (safety grip with red ribbon insert) cat. no. 362 725 ¹⁾ or
- Vacutest Kima cat. no. 149 415 ¹⁾

2.2 Gel sample tubes 100mm

Gel sample tubes are only permitted for the NAT testing of plasma produced from whole blood donation.

When using gel sample tubes, the manufacturer's specifications must be observed in order to avoid haemolytic samples and increased clot formation. Gel sample tubes must be stored in the refrigerator for at least two hours before freezing.

All the sample tubes must be completely filled as the position of the gel layer (among other things) will depend on the filling volume and constitutes a major parameter for the processing.

non ridged, volume 5 ml, with separating gel

Filling height / filling volume: min. 75 mm / corresponds to min. 4.5 ml
max. 85 mm / corresponds to max. 5 ml

- K2E/EDTA K2: Greiner (white safety grip cap with yellow insert) cat. no. 456 058 ¹⁾
- Becton Dickinson (pearly safety grip cap with red insert) cat. no. 362 795 ¹⁾
- Vacutest Kima (white cap) cat. no. 12620 ¹⁾

¹⁾ In case of country-specific deviations of the cat. no. the equivalence of the used sample tube has to be attested by the manufacturer in writing. This confirmation

CONFIDENTIAL: This material is the property of Takeda and must not be disclosed except as authorized in writing.

must be submitted to the contact person at Plasma Logistics & Plasma Analytics on request.


3.0 Barcode position on the sample tube

The beginning of the barcode must be at least 25 mm from the bottom of the tube and the end of the barcode must not be more than 63 mm from the bottom of the tube.


The top-end of the label must not be:

- more than 70 mm from the bottom of the tube (plasma tube 75 mm) (see Picture 1: Position of label/barcode plasma vial).
- more than 80 mm from the bottom of the tube (gel tube 100 mm) (see Picture 2: Position of label/barcode gel vial).

Make sure the labels adhere properly to the tubes and are free of bubbles and/or wrinkles. Make sure the label is straight.

Picture 1: Position of label/ barcode plasma vial

Picture 2: Position of label/ barcode gel vial

4.0 Packing of samples

The NAT sample tubes must be frozen with the caps facing upward and stored at $< -20^{\circ}\text{C}$.

The NAT sample tubes must be put into the labelled box with the nesting insert.

These boxes must be packed into the bigger box. This bigger box must be labelled as well.

Any free space in the bigger box must be filled with filling material.

Either place the bigger box(es) on a separate palette or at the very top of the plasma palette.

No additional paper or electronic documentation must be issued for the delivery of samples that are sent together with plasma shipment.

CONFIDENTIAL: This material is the property of Takeda and must not be disclosed except as authorized in writing.

5.0 Labels for sample tubes and sample boxes

5.1 NAT sample tubes

Length: 30 to 55 mm

Width: 20 to 35 mm

5.1.1 Minimum requirements of content:

- Bleed number in plain text and barcode

5.2 Box for NAT sample tubes

Length: approx. 105 mm

Width: approx. 70 mm

(e.g.: Zweckform# 3426)

5.2.1 Minimum requirements of content:

- Wording: "Samples for NAT testing"
- Supplier name
- Center code
- Shipment number
- Number of samples
- Wording: "Storage and transportation $\leq -20^{\circ}\text{C}$ "

5.3 Bigger box (outer box for NAT sample tube boxes)

Length: approx. 148 mm

Width: approx. 105 mm

(e.g.: Zweckform# 3483)

5.3.1 Minimum requirements of content:

- Wording "Samples for NAT testing"
- Supplier name
- Center code
- Shipment number
- Wording "Storage and transportation $\leq -20^{\circ}\text{C}$ "

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