

ALM-MB-L	Sentrifugehastighet Heraeus Labofuge 400R Prøvetakingspoliklinikken	EQSDocument: 16040 revisjon 1.8 Utskriftsdato: 26.08.2026
Side 1 av 1		

RØRTYPER	Anbefalt G	Koagulere ant. min.	TID ant. min.	Hastighet (RPM)
				Heraeus Labofuge 400R r=16,5 cm
 Gel-rør	1800 – 2200 G	30 min	10	3400
 Li-Heparin gel-rør	2200 G	Romtemp.	10	3400
 EDTA-rør med gel	1800 – 2200 G	Romtemp.	10	3400
 Rør uten gel	1500 – 2200 G	30 min	10	3400
 Urin u/tisetning til mikroskopering	400 – 450 G		5	1500